



Distributed By:



Adelab Scientific
36 Holland Street
Thebarton SA 5031
Ph 08 8234 7955
Fax 08 8234 7897
info@adelab.com.au
www.adelab.com.au

Introducing Accumax Centrifuge Tubes!

Specially designed leak proof screw cap

Large writing area for easy marking & identification

No sharp point or edges


**ISO 9001:2008
ISO 13485:2003**

Easy grip cap for single hand operation even while wearing gloves

Definite graduations for perfect sampling


Unique tamper proof graduation marks resistant to physical & chemical corrosion

Varieties of Accumax Centrifuge Tubes




**50 ML
SELF STANDING**

Withstands upto 7,000 RCF



**50 ML
TUBES**

Withstands upto 20,000 RCF



**15 ML
TUBES**

Withstands upto 17,000 RCF

Specifications

- Manufactured from FDA approved resin
- Certified RNase, DNase, Human DNA, Endotoxin & Pyrogen free
- Free from heavy metals & additives
- Manufactured in a fully automatic, human touch free plant with class 100K clean room facility
- Can be centrifuged at up to 20,000xg RCF
- Withstand temperatures from - 80 to +121degC
- Fit all standard and high capacity centrifuge rotors
- Supplied in double safety self-standing bulk bags & Eco friendly Cardboard Racks

15 ML & 50 ML Centrifuge Tubes in Eco Friendly Rack



15 ML & 50 ML Individually Wrapped Centrifuge Tubes



Ordering Information

CAT NO.	Description	Qty	Unit	Price Ex-GST
ACT-15-B-S	Accumax 15 ml Centrifuge Tube with flat Cap, Sterile, Max RCF:17,000xg, pkt/500	1	pkt	\$58.00
ACT-50-B	Accumax 50 ml Centrifuge Tube with flat Cap, Max RCF: 20,000xg, pkt/500	1	pkt	\$63.00
ACT-50-B-S	Accumax 50 ml Centrifuge Tube with flat Cap, Sterile, Max RCF: 20,000xg, pkt/500	1	pkt	\$73.00
ACT-50-K-B-S	Accumax 50 ml Self Standing Centrifuge Tube with flat Cap, Sterile, Max RCF: 7,000xg, pkt/500	1	pkt	\$80.00