

Adelab's pH Meters

Bench

\$872.00
Ex-GST



EUTECH-pH 700 Bench Meter

Includes pH Electrode (ECFC7252101B), ATC Probe (ECPHTEMPB01P), Stand, AC Adapter.



CP-511 pH, Temp

\$829.00
Ex-GST



- Accurate pH, mV (redox potential) and temp measurements.
- Calibration of the pH electrode: 1-3 points.
- Automatic buffer detection 4.00; 7.00; 10.00 pH.
- Automatic temperature compensation.



\$675.00
Ex-GST



Starter 3100

\$405.00
Ex-GST



Starter 2100

OHAUS Starter 2100

Starter 2100 is the bench meter of choice for those looking for a cost-effective pH meter that will provide reliable measurements time and time again.

OHAUS Starter 3100

Starter 3100 offers many premier features that propel the bench meter to an advanced level of operational efficiency.

OHAUS Starter 5000

Starter 5000 provides accurate pH and ORP measurement for high-level experiments and research. Enhanced features, advanced technology, and high performance have been combined to support complex laboratory tasks requiring pH measurement.

\$1,170.00
Ex-GST



Starter 5000

Distributed By:



Adelab Scientific
36 Holland Street
Thebarton SA 5031
Ph 08 8234 7955
Fax 08 8234 7897
info@adelab.com.au
www.adelab.com.au

Portable Pocket PH Meters

Ionode-CPX-105

- Includes Ionode's famous IJ-44 pH Sensor
- the right choice for difficult samples
- Suitable in food, dairy, meat, waste-water and solids
- toughened spear membrane for direct penetration
- double Junction with removable sleeve defeats contamination
- up to 200 hours of continual use from battery power
- 1 to 3 point calibration
- automatic detection of buffer solutions

\$108.00
Ex-GST



Eutech EcoTestr pH 2 Designed for fuss-free pH measurements on the go. Accurate and easy to use.



\$386.00
Ex-GST



\$205.00
Ex-GST

Eutech pHTestr® 30

Uniquely designed to last longer with more user-friendly features and lower operating cost requirement. Double junction electrode makes the testers ideal for dirty water and solutions with heavy metals or contaminants, while increased reference gel volume offers significantly longer electrode life.

PH Electrodes with BNC connector



Eutech-ECFC7252101B: General Purpose Plastic-body Single Junction Gel-Filled pH Combination Electrode, Range; 0-13 PH, 12x90mm, BNC connector, 1m cable-**Price \$140.00 Ex-GST.**



Sensorex-S200C/BNC Epoxy Plastic Body General Purpose pH Electrode with BNC Connector, 0 to 14 pH, 12mm Diameter x 150mm Length-**Price \$157.00 Ex-GST.**



IONODE-IJ44C Intermediate Junction pH Electrode, pH 0 - 14 with replaceable electrolyte, general purpose, suitable for insertion into meat, cheese -**Price \$185.00 Ex-GST.**

Adelab's PH- Buffers, Storage and Cleaning solutions.

SCHS2004	Scharlau Buffer pH 4, red coloured 500ml, NIST traceable, certified ± 0.01 pH units-each	\$16.00 Ex-Gst
SCHS2007	Scharlau Buffer pH 7, yellow coloured 500ml, NIST traceable, certified ± 0.01 pH units-each	\$16.00 Ex-Gst
SCHS2010	Scharlau Buffer pH 10, blue coloured 500ml, NIST traceable, certified ± 0.01 pH units-each	\$16.00 Ex-Gst
AEEC-AS-ST-S1-05	Sitest Storage Solution for pH electrodes (500mL)-each	\$50.00 Ex-Gst
AEEC-DPC-BT	Eutech 480 mL Protein Cleaning Solution for pH electrode (480mL)-each	\$45.00 Ex-Gst