

Specials

THERMOMETERS AND TIMER

(While stocks last)

Adelab
Scientific



DIGITAL TIMER

12 or 24 hour clock

Count up or Count down in hours, minutes and seconds

\$ 15.00 (excl. GST)

DIGITAL THERMOMETER

Range: -50 to +70°C

-58 to 158°F


Resolution: 0.1°C

Accuracy: 0.1°C

Cable: 1 metre

Battery: AAA (incl.)

\$ 15.00 (excl. GST)



DIGITAL THERMOMETER with MAX-MIN function

Range: -50 to +200°C

-58 to +392°F

Resolution: 0.1°C

Probe: 125mm

Battery: LR44(incl.)

\$ 15.00 (excl. GST)



DIGITAL THERMOMETER with MAX-MIN function

Range: -50 to +200°C

-58 to +392°F

Resolution: 0.1°C

Probe: 125mm

Battery: A76 (incl.)

\$ 15.00 (excl. GST)



Specials

MAGNETIC STIRRING BARS, PTFE COATED

**COWIE
BRAND**

Guaranteed inert and leak-free
Isostatically molded in pure, FDA approved, PTFE
Individually polished to stay cleaner for longer

| | | | Qty | Each |
|--------------|---------|-----------------------------------|-----|--------|
| LT-001.106 | Cowie | Stirring Bar Cylindrical 10 x3mm | 87 | \$2.50 |
| LT-001.112 | Cowie | Stirring Bar Cylindrical 12x4.5mm | 471 | \$2.50 |
| LT-001.113.6 | Cowie | Stirring Bar Cylindrical 13x6mm | 27 | \$2.50 |
| LT-001.115 | Cowie | Stirring Bar Cylindrical 15x4.5mm | 63 | \$2.50 |
| LT-001.120 | Cowie | Stirring Bar Cylindrical 20x6mm | 11 | \$2.50 |
| LT-001.125 | Cowie | Stirring Bar Cylindrical 25x6mm | 52 | \$2.50 |
| LT-001.130 | Cowie | Stirring Bar Cylindrical 30x6mm | 14 | \$2.50 |
| KAR758 | Kartell | Stirring bar Cylindrical 30x6mm | 500 | \$2.50 |
| LT-001.135 | Cowie | Stirring Bar Cylindrical 35x6mm | 126 | \$2.50 |
| LT-001.160 | Cowie | Stirring Bar Cylindrical 60x10mm | 29 | \$2.50 |



| | | | | |
|--------------|-------|-----------------------------|----|--------|
| LT-001.312 | Cowie | Stirring Bar Pivot 12x4.5mm | 60 | \$2.50 |
| LT-001.315.8 | Cowie | Stirring Bar Pivot 15x8mm | 11 | \$2.50 |
| LT-001.320 | Cowie | Stirring Bar Pivot 20x6mm | 48 | \$2.50 |
| LT-001.330 | Cowie | Stirring Bar Pivot 30x6mm | 38 | \$2.50 |
| LT-001.340 | Cowie | Stirring Bar Pivot 40x8mm | 2 | \$2.50 |
| LT-001.345 | Cowie | Stirring Bar Pivot 45x5mm | 4 | \$2.50 |
| LT-001.350 | Cowie | Stirring Bar Pivot 50x8mm | 14 | \$2.50 |
| LT-001.360 | Cowie | Stirring Bar Pivot 60x10mm | 39 | \$2.50 |
| LT-001.370 | Cowie | Stirring Bar Pivot 70x10mm | 6 | \$2.50 |



| | | | | |
|------------|-------|--------------------------------|----|--------|
| LT-001.513 | Cowie | Stirring Bar Octahedral 13x8mm | 10 | \$2.50 |
|------------|-------|--------------------------------|----|--------|



| | | | | |
|------------|-------|--------------------------------|-----|--------|
| LT-001.610 | Cowie | Stirring Bar Egg Shape 10x5mm | 129 | \$2.50 |
| LT-001.615 | Cowie | Stirring Bar Egg Shape 15x6mm | 464 | \$2.50 |
| LT-001.620 | Cowie | Stirring Bar Egg Shape 20x10mm | 20 | \$2.50 |
| LT-001.625 | Cowie | Stirring Bar Egg Shape 25x12mm | 24 | \$2.50 |
| LT-001.630 | Cowie | Stirring Bar Egg Shape 30x16mm | 9 | \$2.50 |
| LT-001.635 | Cowie | Stirring Bar Egg Shape 35x16mm | 17 | \$2.50 |
| LT-001.640 | Cowie | Stirring Bar Egg Shape 40x20mm | 13 | \$2.50 |



| | | | | |
|-------------|-------|--|----|---------|
| LT-003.300 | Cowie | Stirring Rod 300mm x 6mm PTFE coated steel | 17 | \$8.00 |
| LT-004.1350 | Cowie | Stirring Bar Retriever polypropylene 350mm | 40 | \$10.00 |



ONLY WHILE STOCKS LAST