



**ADELAB SCIENTIFIC SUPPLIES AND  
SUPPORTS THE WATER  
PURIFICATION GROUP OF MILLIPORE'S  
LIFE SCIENCES DIVISION**

**ENCLOSED IS A SUMMARY OF PRODUCTS  
AVAILABLE TO SUIT A WIDE RANGE OF  
PURIFIED WATER APPLICATIONS.  
IF YOU SHOULD WISH TO OBTAIN  
FURTHER INFORMATION ON ANY OF  
THESE, PLEASE ASK ANY OF OUR STAFF**

**WE THANK YOU FOR YOUR INTEREST**

---

## Water Types & Applications

---

**Type III:** glassware rinsing, filling of heating baths and autoclaves, recirculating cooling systems as well as a feed supply to Type I Laboratory Water Systems



Water baths & circulators



Rotary evaporators



Autoclaves

**Type II:** buffers, pH solutions, microbiological culture media preparation, feed for clinical analysers and weatherometers, preparation of reagents for chemical analysis or synthesis



pH buffers



Media preparation




Titrators



Glassware washers

**Type I:** HPLC mobile phase preparation, blanks and sample dilution in GC, HPLC, AA, ICP-MS and other advanced analytical techniques, preparation of buffers and culture media for mammalian cell culture, IVF, reagents for molecular biology and other critical requirements



Molecular biology



Analytical instruments



HPLC



Cell culture

Type 3 water - for basic laboratory applications

Model RiOs 3  
Producing 3 L/hr



Models

RiOs 5	5 L/hr
RiOs 8	8 L/hr
RiOs 16	16 L/hr

Models

RiOs 30	30 L/hr
RiOs 50	50 L/hr
RiOs 100	100 L/hr
RiOs 150	150 L/hr
RiOs 200	200 L/hr




Type 2 water - for general laboratory applications



Models

- ELIX 3 ADVANTAGE 3 L/hr
- ELIX 5 ADVANTAGE 5 L/hr
- ELIX 10 ADVANTAGE 10L/hr
- ELIX 15 ADVANTAGE 15 L/hr
- ELIX 3 3 L/hr
- ELIX 5 5 L/hr
- ELIX 10 10L/hr
- ELIX 15 15 L/hr

All ELIX models utilise Millipore's advanced Electro-deionization technology



Models

- ELIX 20 20 L/hr
- ELIX 35 35 L/hr
- ELIX 70 70 L/hr
- ELIX 100 100 L/hr



Purified Water Storage


RESERVOIRS FOR STORING PURIFIED TYPES 3 AND 2 WATER SEALED FROM ATMOSPHERIC AND ENVIRONMENTAL CONTAMINATION



Reservoirs  
30, 60 & 100 Litres



Auto UV sanitization module



Vent Filters

Reservoirs  
250 & 350 Litres



Type 1 water - for critical laboratory applications



**Milli-Q Simplicity & Milli-Q Simplicity UV**



**Milli-Q Synergy & Milli-Q Synergy UV**



**Milli-Q Advantage A10 with true TOC monitor with up to 3 Q-Pods**

**Milli-Q Reference and Milli-Q Reference+ (with TOC indication)**



Models require a supply of pre-purified water

Type 1 water - for critical laboratory applications

Ultra-pure water direct from your potable tap water supply !



Milli-Q Direct Q 3

Milli-Q Direct Q 8  
Milli-Q Direct Q 16



Milli-Q Integral System



Models available:

Integral 3, Integral 5, Integral 10, Integral 15

All Millipore Laboratory Water Division products are available in South Australia from

## ADELAB SCIENTIFIC

36 Holland Street, Thebarton SA 5031

Telephone: 8234 7955 Fax: 8234 7897

Email: [info@adelab.com.au](mailto:info@adelab.com.au) Web: [www.adelab.com.au](http://www.adelab.com.au)



### SYSTEMS

### SERVICE



### CONSUMABLES

Call Adelab Scientific to order replacement cartridges, arrange a service visit or discuss your laboratory's water requirements.