


# FavorPrep™ Stool DNA Isolation Mini Kit



## Features

- **Time Saving:** Rapid isolation of ready-to-use DNA within 60 minutes without phenol/chloroform extraction.
- **High purity:** Eliminate humic acid, polysaccharides, phenol compounds, and enzyme inhibitor from stool sample.
- **Sample Size:** 50~200mg of fresh or frozen stool sample.
- **Format:** Spin Column



## Specifications

**Sample:** 50~200mg of fresh or frozen stool sample.

**Binding Capacity:** Up to 30µg genomic DNA

**Operation Time:** Within 60 minutes

## Order Information

Cat. No.	Product Name	Size
<b>FASTI 001</b>	FavorPrep™ Stool DNA Isolation	50 preps
<b>FASTI 001-1</b>	Mini Kit	100 preps

### Kit Component

Glass Beads SDE1 Buffer SDE2 Buffer SDE3 Buffer  
SDEX Buffer Wash Buffer (concentrated)  
Elution Buffer Proteinase K SDE Mini Column  
2.0ml Tube (for Collection)